

Product Datasheet

Human S100A12 ELISA Kit (orb1809166)

| | |
|------------------------|--|
| Description | Human S100A12 ELISA Kit |
| Reactivity | Human |
| Range | 78.1-5000pg/mL |
| Storage | Store at 4°C & -20°C for 12 months. |
| Note | For research use only |
| Sample Types | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Sensitivity | 34pg/mL |
| Expiration Date | Please enquire. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)