


## Product Datasheet

### ARN509 (orb180925)

<b>Description</b>	ARN-509 is a selective and competitive androgen receptor inhibitor with IC50 of 16 nM, useful for pr...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{21}H_{15}F_4N_5O_2S$
<b>MW</b>	477.4347
<b>CAS Number</b>	[956104-40-8]
<b>SMILES</b>	<chem>S=C1N(C2C([H])=NC(C#N)=C(C(F)(F)F)C=2[H])C(C2(C([H]))([H])C([H]))([H])C2([H])[H])N1C1C([H])=C([H])C(C(N([H])C([H]))([H])[H])=O)=C(C=1[H])F)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.



Chemical structure of ARN509