

## Product Datasheet

### Human AngII ELISA Kit (orb1809286)

<b>Description</b>	Human AngII ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	3.12-200pg/mL
<b>Storage</b>	Store at 4°C & -20°C for 12 months.
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma or other biological fluids.
<b>Sensitivity</b>	1.13pg/mL
<b>Expiration Date</b>	Please enquire.