

## Product Datasheet

### Mouse a2M ELISA Kit (orb1809304)

<b>Description</b>	Mouse a2M ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	9.37-600ng/mL
<b>Storage</b>	Store at 4°C & -20°C for 12 months.
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, urine or other biological fluids.
<b>Sensitivity</b>	3.7ng/mL
<b>Expiration Date</b>	Please enquire.