

## Product Datasheet

### Human SPD ELISA Kit (orb1809322)

<b>Description</b>	Human SPD ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	6.25-400ng/mL
<b>Storage</b>	Store at 4°C & -20°C for 12 months.
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates, cell lysates, lung lavage fluid, cell culture supernates or other biological fluids.
<b>Sensitivity</b>	2.16ng/mL
<b>Expiration Date</b>	Please enquire.