

## Product Datasheet

### Mouse ANGPTL3 ELISA Kit (orb1809385)

<b>Description</b>	Mouse ANGPTL3 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.156-10ng/mL
<b>Storage</b>	Store at 4°C & -20°C for 12 months.
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates or other biological fluids.
<b>Sensitivity</b>	0.061ng/mL
<b>Expiration Date</b>	Please enquire.