


Product Datasheet

Bi 2536 (orb180942)

Description	BI2536 is a potent PI3K inhibitor with IC50 of 0.83 nM, BI2536 shows 4- and 11-fold greater selectiv...
Note	For research use only
Formula	C ₂₈ H ₃₉ N ₇ O ₃
MW	521.6544
CAS Number	[755038-02-9]
SMILES	<chem>O=C1[C@@]([H])(C([H])([H])C([H])([H])[H])N(C2C(=C([H])N=C(N([H])C3C([H])=C([H])C(=C([H])C=3OC([H])([H])[H])C(N([H])C3([H])C([H])([H])C([H])([H])N(C([H])([H])[H])C([H])([H])C3([H])[H])=O)N=2)N1C([H])([H])[H])C1([H])C([H])([H])C([H])([H])C([H])([H])C1([H])[</chem>
Expiration Date	12 months from date of receipt.



Chemical structure of Bi 2536