

Product Datasheet

rabbit IgG Antibody (BF594) (orb1809473)

| | |
|----------------------------|---|
| Description | rabbit IgG Antibody (BF594) |
| Reactivity | Rabbit |
| Conjugation | BF594 |
| Tested Applications | IF |
| Target | rabbit IgG |
| Preservatives | 10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol. |
| Concentration | 1mg/ml |
| Storage | Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |