

## Product Datasheet

### SARS-CoV-2 Antibody (HRP) (orb1809474)

<b>Description</b>	SARS-CoV-2 Antibody (HRP)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Recombinant SARS-CoV-2 Nucleocapsid protein
<b>Target</b>	SARS-CoV-2
<b>Preservatives</b>	0.2 µm filtered solution of PBS, pH 7.4.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1/kappa
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P0DTC2</a>
<b>Dilution Range</b>	WB=1:500-2000
<b>Expiration Date</b>	12 months from date of receipt.