


## Product Datasheet

### BX912 (orb180950)

<b>Description</b>	BX-912 is a selective inhibitor of 3-Phosphoinositide-dependent Kinase-1 (PDK1) with IC50 of 12 nM; ...
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>20</sub> H <sub>23</sub> BrN <sub>8</sub> O
<b>MW</b>	471.35
<b>CAS Number</b>	[702674-56-4]
<b>SMILES</b>	BrC1=C(NCCC2=CNC=N2)N=C(NC3=CC(NC(N4CCCC4)=O)=CC=C3)N=C1
<b>Expiration Date</b>	12 months from date of receipt.



Chemical structure of BX912