

Product Datasheet

SARS-CoV-2 Antibody (HRP) (orb1809539)

Description	SARS-CoV-2 Antibody (HRP)
Species/Host	Mouse
Reactivity	Virus
Conjugation	HRP
Tested Applications	WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protein (1-419/419aa)
Target	SARS-CoV-2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	PODTC2
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.