

Product Datasheet

SARS-CoV-2 Antibody (Biotin) (orb1809689)

| | |
|----------------------------|--|
| Description | SARS-CoV-2 Antibody (Biotin) |
| Species/Host | Mouse |
| Reactivity | Virus |
| Conjugation | Biotin |
| Tested Applications | WB |
| Immunogen | Recombinant SARS-CoV-2 Nucleocapsid Protein |
| Target | SARS-CoV-2 |
| Preservatives | PBS, pH 7.4 |
| Concentration | 3.5mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Monoclonal |
| Uniprot ID | P0DTC2 |
| Expiration Date | 12 months from date of receipt. |