

Product Datasheet

Talnetant (SB 223412) (orb180970)

Description	Talnetant (SB 223412) is a potent and selective NK3 receptor antagonist($k_i=1.4$ nM, hNK-3-CHO); 100-f...
Note	For research use only
Formula	$C_{25}H_{22}N_2O_2$
MW	382.45438
CAS Number	[174636-32-9]
SMILES	<chem>O=C(C1=C(O)C(C2=CC=CC=C2)=NC3=CC=CC=C13)N[C@H](C4=CC=CC=C4)CC</chem>
Expiration Date	12 months from date of receipt.