

## Product Datasheet

### PKI402 (orb180983)

<b>Description</b>	PKI-402 is a potent dual pan-PI3K/mTOR inhibitor targeting PI3K $\alpha/\beta/\gamma/\delta$ and mTOR with IC50 of 2 n...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{29}H_{34}N_{10}O_3$
<b>MW</b>	570.65
<b>CAS Number</b>	[1173204-81-3]
<b>SMILES</b>	<chem>CCN1N=NC2=C1N=C(N=C2N3CCOCC3)C4=CC=C(C=C4)NC(NC5=CC=C(C=C5)C(N6CCN(CC6)C)=O)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.