

## Product Datasheet

### Dacinostat (NVP-LAQ824, LAQ824) (orb180988)

**Description** Dacinostat (NVP-LAQ824, LAQ824) is a highly potent HDAC inhibitor with IC50 of 32 nM.

**Note** For research use only


**Formula** C<sub>22</sub>H<sub>25</sub>N<sub>3</sub>O<sub>3</sub>

**MW** 379.4522

**CAS Number** [404951-53-7]

**SMILES** O([H])C([H])([H])C([H])([H])N(C([H])([H])C1C([H])=C([H])C(/C([H])=C(\[H])/C(N([H])O[H])=O)=C([H])C=1[H])C([H])([H])C([H])([H])C1=C([H])N([H])C2=C([H])C([H])=C([H])C([H])=C12

**Expiration Date** 12 months from date of receipt.



Chemical structure of Dacinostat (NVP-LA...