

Product Datasheet

Tarafenacin(SVT-40776) (orb180991)

Description	Tarafenacin(SVT-40776) is a highly selective M3 muscarinic receptor antagonist (Ki= 0.19 nM), ~200 f...
Note	For research use only
Formula	$C_{25}H_{26}F_4N_2O_8$
MW	558.476153
CAS Number	[385367-47-5]
SMILES	<chem>O=C(O[C@H]1CN2CCC1CC2)N(C3=CC=CC(F)=C3)CC4=CC(F)=C(F)C(F)=C4</chem>
Expiration Date	12 months from date of receipt.