

Product Datasheet

IKBKB Antibody (BF750) (orb1809935)

| | |
|----------------------------|---|
| Description | IKBKB Antibody (BF750) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | BF750 |
| Tested Applications | ICC, IF |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from human IKK beta around the phosphorylation site of S474 MA(p-S)MS |
| Target | IKBKB |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | O14920 |
| Expiration Date | 12 months from date of receipt. |