

Product Datasheet

IL2RA Antibody (PE/Cy5.5) (orb1810079)

Description	IL2RA Antibody (PE/Cy5.5)
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy5.5
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthetic peptide derived from human CD25
Target	IL2RA
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Concentration	1mg/ml
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG2a, k
Clonality	Monoclonal
Uniprot ID	P01590
Expiration Date	12 months from date of receipt.