


Product Datasheet

IC-87114 (orb181021)

Description	IC-87114 was the first isoform-selective PI3K inhibitor. p110 β (IC ₅₀ = 0.13 μ M) vs. p110 α (IC ₅₀ = 2...
Note	For research use only
Formula	C ₂₂ H ₁₉ N ₇ O
MW	397.43256
CAS Number	[371242-69-2]
SMILES	O=C1N(C(CN2C3=C(C(N)=NC=N3)N=C2)=NC4=CC=CC(C)=C14)C5=C(C)C=CC=C5
Expiration Date	12 months from date of receipt.



Chemical structure of IC-87114