


Product Datasheet

LGK-974 (orb181048)

Description	LGK-974 is a potent and specific PORCN inhibitor, and inhibits Wnt signaling with IC50 of 0.4 nM.
Note	For research use only
Formula	C ₂₃ H ₂₀ N ₆ O
MW	396.44
CAS Number	[1243244-14-5]
SMILES	O=C(NC1=NC=C(C2=NC=CN=C2)C=C1)CC3=CN=C(C4=CC(C)=NC=C4)C(C)=C3
Expiration Date	12 months from date of receipt.



Chemical structure of LGK-974