

Product Datasheet

RB1 Antibody (BF750) (orb1810671)

Description	RB1 Antibody (BF750)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF750
Tested Applications	ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Rb around the phosphorylation site of Thr821 LP(p-T)PT
Target	RB1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P06400
Expiration Date	12 months from date of receipt.