

Product Datasheet

Monkeypox virus Antibody (HRP) (orb1810705)

| | |
|----------------------------|---|
| Description | Monkeypox virus Antibody (HRP) |
| Species/Host | Rabbit |
| Reactivity | Virus |
| Conjugation | HRP |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant MPXV M1R protein (1-182/250aa) |
| Target | Monkeypox virus |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q8V502 |
| Expiration Date | 12 months from date of receipt. |