

Product Datasheet

Monkeypox virus Antibody (Biotin) (orb1810706)

Description	Monkeypox virus Antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant MPXV M1R protein (1-182/250aa)
Target	Monkeypox virus
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8V502
Expiration Date	12 months from date of receipt.