

## Product Datasheet

### OSI 930 (orb181074)

|                        |   |
|------------------------|---|
| <b>Description</b>     | OSI-930 is a potent inhibitor of Kit, KDR and CSF-1R with IC50 of 80 nM, 9 nM and 15 nM, respectivel... |
| <b>Note</b>            | For research use only   |
| <b>Formula</b>         | $C_{22}H_{16}F_3N_3O_2S$  |
| <b>MW</b>              | 443.44  |
| <b>CAS Number</b>      | [728033-96-3]   |
| <b>SMILES</b>          | <chem>O=C(C1=C(NCC2=CC=NC3=CC=CC=C23)C=CS1)NC4=CC=C(OC(F)(F)F)C=C4</chem>                               |
| <b>Expiration Date</b> | 12 months from date of receipt.   |