

Product Datasheet

HPV16 Antibody (APC/Cy5.5) (orb1810745)

Description	HPV16 Antibody (APC/Cy5.5)
Species/Host	Rabbit
Reactivity	Human
Conjugation	APC/Cy5.5
Tested Applications	FC
Immunogen	Recombinant HPV16 E6
Target	HPV16
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03129
Expiration Date	12 months from date of receipt.