

Product Datasheet

PI103 (orb181087)

Description	PI103 is a potent inhibitor with low IC50 values against recombinant PI3K isoforms p110alpha (IC50= ...
Note	For research use only
Formula	$C_{19}H_{16}N_4O_3$
MW	348.36
CAS Number	[371935-74-9]
SMILES	<chem>OC1=CC(C2=NC3=C(C(N4CCOCC4)=N2)OC5=C3C=CC=N5)=CC=C1</chem>
Expiration Date	12 months from date of receipt.