

## Product Datasheet

### PIK-294 (orb181089)

<b>Description</b>	PIK-294 is a highly selective p110δ inhibitor with IC50 of 10 nM, 1000-, 49- and 16-fold less poten...
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>28</sub> H <sub>23</sub> N <sub>7</sub> O <sub>2</sub>
<b>MW</b>	489.53
<b>CAS Number</b>	[900185-02-6]
<b>SMILES</b>	<chem>N1C2=C(C(C)=CC=C2)C(=O)N(C2=CC=CC=C2)C=1CN1C2C(C(C3=CC=CC(O)=C3)=N1)=C(N)N=CN=2</chem>
<b>Expiration Date</b>	12 months from date of receipt.