

## Product Datasheet

### PIK-75 (orb181090)

<b>Description</b>	PIK-75 is a p110 $\alpha$ inhibitor with IC <sub>50</sub> of 5.8 nM (200-fold more potently than p110 $\beta$ ), isoform-speci...
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>16</sub> H <sub>15</sub> BrClN <sub>5</sub> O <sub>4</sub> S
<b>MW</b>	488.741
<b>CAS Number</b>	[372196-77-5]
<b>SMILES</b>	CC1=C(C=C(C=C1)[N+](=O)[O-])S(=O)(=O)N(C)/N=C/C2=CN=C3N2C=C(C=C3)Br.Cl
<b>Expiration Date</b>	12 months from date of receipt.