

Product Datasheet

Gedatolisib(PF-05212384) (orb181093)

Description	PKI-587 is a highly potent dual PI3K/mTOR kinase inhibitor with IC50 of 0.4 nM and < 0.1 μM for PI3...
Note	For research use only
Formula	C ₃₂ H ₄₁ N ₉ O ₄
MW	615.73
CAS Number	[1197160-78-3]
SMILES	CN(C)C1CCN(CC1)C(=O)C2=CC=C(C=C2)NC(=O)NC3=CC=C(C=C3)C4=NC(=NC(=N4)N5CCOCC5)N6CCOCC6
Expiration Date	12 months from date of receipt.