

Product Datasheet

PF-8380 (orb181098)

Description	PF-8380 is a potent and specific autotaxin inhibitor with an IC50 value of 2.8 nM.
Note	For research use only
Formula	$C_{22}H_{21}Cl_2N_3O_5$
MW	478.33
CAS Number	[1144035-53-9]
SMILES	<chem>C1CN(CCN1CCC(=O)C2=CC3=C(C=C2)NC(=O)O3)C(=O)OCC4=CC(=CC(=C4)Cl)Cl</chem>
Expiration Date	12 months from date of receipt.