

Product Datasheet

B3GAT1 Antibody (PE/Cy5.5) (orb1811128)

Description	B3GAT1 Antibody (PE/Cy5.5)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human B3GAT1 (221-320/334aa)
Target	B3GAT1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9P2W7
Dilution Range	Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.