

## Product Datasheet

### SF1670 (orb181129)

<b>Description</b>	SF1670 is a selective PTEN (phosphatase and tensin homolog) inhibitor with an IC50 of 2 micromolar.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{19}H_{17}NO_3$
<b>MW</b>	307.34
<b>CAS Number</b>	[345630-40-2]
<b>SMILES</b>	<chem>CC(C)(C)C(=O)NC1=CC2=C(C=C1)C3=CC=CC=C3C(=O)C2=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.