

## Product Datasheet

### Pracinostat(SB939) (orb181135)

<b>Description</b>	Pracinostat(SB939) is a potent pan-HDAC inhibitor with IC50 of 40-140 nM with exception for HDAC6; h...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{20}H_{30}N_4O_2$
<b>MW</b>	358.48
<b>CAS Number</b>	[929016-96-6]
<b>SMILES</b>	<chem>C(NO)(=O)/C=C/C1=CC=C2C(=C1)N=C(CCCC)N2CCN(CC)CC</chem>
<b>Expiration Date</b>	12 months from date of receipt.