

## Product Datasheet

### TAK 715 (orb181153)

<b>Description</b>	TAK-715 is a p38 MAPK inhibitor for p38 $\alpha$ with IC50 of 7.1 nM, 28-fold more selective for p38 $\alpha$ over...
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>24</sub> H <sub>21</sub> N <sub>3</sub> O <sub>5</sub>
<b>MW</b>	399.51
<b>CAS Number</b>	[303162-79-0]
<b>SMILES</b>	<chem>O=C(C1=CC=CC=C1)NC2=NC=CC(C3=C(N=C(S3)CC)C4=CC=CC(C)=C4)=C2</chem>
<b>Expiration Date</b>	12 months from date of receipt.