


## Product Datasheet

### Tenovin-3 (orb181160)

<b>Description</b>	Tenovin-3 is a small molecule activator of p53 transcriptional activity.
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>18</sub> H <sub>21</sub> N <sub>3</sub> OS
<b>MW</b>	327.44384
<b>CAS Number</b>	[1011301-27-1]
<b>SMILES</b>	<chem>CC(C)(C)C1=CC=C(C=C1)C(=O)NC(=S)NC2=CC=C(C=C2)N</chem>
<b>Expiration Date</b>	12 months from date of receipt.



Chemical structure of Tenovin-3