

## Product Datasheet

### Tenovin-6 (orb181161)

<b>Description</b>	Tenovin-6 is the water soluble analog of Tenovin-1 and acts as a potent SIRT1 (IC50=21 uM) and SIRT2...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{25}H_{34}N_4O_2S$
<b>MW</b>	454.6
<b>CAS Number</b>	[1011557-82-6]
<b>SMILES</b>	<chem>CC(C)(C)C1=CC=C(C=C1)C(=O)NC(=S)NC2=CC=C(C=C2)NC(=O)CCCCN(C)C</chem>
<b>Expiration Date</b>	12 months from date of receipt.