


Product Datasheet

3-Methyladenine (orb181186)

Description	3-Methyladenine (3-MA) is a selective PI3K inhibitor for Vps34 and PI3Kγ with IC50 of 25 μM and 60...
Note	For research use only
Formula	C ₆ H ₇ N ₅
MW	149.1533
CAS Number	[5142-23-4]
SMILES	<chem>N1(C([H])([H])[H])C([H])=N/C(/C2=C1N=C([H])N2[H])=N[H]</chem>
Expiration Date	12 months from date of receipt.



Chemical structure of 3-Methyladenine