

Product Datasheet

CD274 Antibody (PE/Cy7) (orb1812076)

Description	CD274 Antibody (PE/Cy7)
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy7
Tested Applications	FC
Target	CD274
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Isotype	Mouse IgG2b, k
Clonality	Monoclonal
Uniprot ID	Q9NZQ7
Expiration Date	12 months from date of receipt.