

Product Datasheet

HS-173 (orb181255)

Description	HS-173 is a potent PI3K α inhibitor with IC ₅₀ of 0.8 nM.
Note	For research use only
Formula	C ₂₁ H ₁₈ N ₄ O ₄ S
MW	422.46
CAS Number	[1276110-06-5]
SMILES	O=C(C1=CN=C2C=CC(C3=CC(NS(=O)(C4=CC=CC=C4)=O)=CN=C3)=CN21)OCC
Expiration Date	12 months from date of receipt.