

Product Datasheet

Nano Tag Antibody (Biotin) (orb1812605)

Description	Nano Tag Antibody (Biotin)
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide Nano-Tag (9)
Target	Nano Tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.