

Product Datasheet

IWP-L6 (orb181262)

Description	IWP L6 is a Porcn inhibitor with EC50 of 0.5 nM.
Note	For research use only
Formula	$C_{25}H_{20}N_4O_2S_2$
MW	472.58
CAS Number	[1427782-89-5]
SMILES	<chem>O=C(NC1=NC=C(C2=CC=CC=C2)C=C1)CSC3=NC(CCS4)=C4C(N3C5=CC=CC=C5)=O</chem>
Expiration Date	12 months from date of receipt.