

Product Datasheet

TUBG1 Antibody (Biotin) (orb1812631)

Description	TUBG1 Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Gamma Tubulin
Target	TUBG1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P23258
Expiration Date	12 months from date of receipt.