

Product Datasheet

rbcl Antibody (Biotin) (orb1812692)

Description	rbcl Antibody (Biotin)
Species/Host	Mouse
Reactivity	A. thaliana
Conjugation	Biotin
Tested Applications	WB
Immunogen	D-Ribulose 1,5-Diphosphate Carboxylase purified
Target	rbcl
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P00875
Expiration Date	12 months from date of receipt.