

Product Datasheet

KU-60019 (orb181274)

Description	KU-60019 is an improved analogue of KU-55933, with IC50 of 6.3 nM for ATM.
Note	For research use only
Formula	$C_{30}N_3O_5SH_{33}$
MW	547.6651
CAS Number	[925701-46-8]
SMILES	<chem>O=C(NC1=CC(CC2=C(C(C3=CC(C=C(N4CCOCC4)O3)=O)=CC=C2)S5)=C5C=C1)CN6C[C@H](C)O[C@H](C)C6</chem>
Expiration Date	12 months from date of receipt.