

Product Datasheet

ML-323 (orb181292)

Description	ML323 is a reversible, potent USP1-UAF1 inhibitor with IC50 of 76 nM in a Ub-Rho assay and 174 nM an...
Note	For research use only
Formula	$C_{23}H_{24}N_6$
MW	384.48
CAS Number	[1572414-83-5]
SMILES	<chem>C(C1=CC=CC=C1C(C)C)1=NC=C(C)C(NCC2=CC=C(N3C=CN=N3)C=C2)=N1</chem>
Expiration Date	12 months from date of receipt.