

Product Datasheet

TG100-115 (orb181334)

Description	TG100115 is a PI3Ky/6 inhibitor with IC50 of 83 nM/235 nM, with little effect on PI3K α / β . Phase ...
Note	For research use only
Formula	C ₁₈ H ₁₄ N ₆ O ₂
MW	346.34
CAS Number	[677297-51-7]
SMILES	NC1=C(C(N=C2C3=CC(O)=CC=C3)=NC(N)=N1)N=C2C4=CC(O)=CC=C4
Expiration Date	12 months from date of receipt.