

Product Datasheet

SARS-CoV-2 Antibody (Biotin) (orb1813350)

Description	SARS-CoV-2 Antibody (Biotin)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 (2019-nCoV) Nucleocapsid protein full length
Target	SARS-CoV-2
Preservatives	0.01M PBS(pH7.4) with 0.03% Proclin300
Concentration	≥1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P0DTC2
Expiration Date	12 months from date of receipt.