

Product Datasheet

HPV16 Antibody (BF647) (orb1813425)

| | |
|----------------------------|--|
| Description | HPV16 Antibody (BF647) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | BF647 |
| Tested Applications | FC, IF |
| Immunogen | KLH conjugated synthetic peptide derived from human HPV16 E6 (85-158/158aa) |
| Target | HPV16 |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P03129 |
| Expiration Date | 12 months from date of receipt. |