

Product Datasheet

HPV16 Antibody (PerCP/Cy5.5) (orb1813461)

Description	HPV16 Antibody (PerCP/Cy5.5)
Species/Host	Rabbit
Reactivity	Human, Virus
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human HPV16 E6+ HPV18 E6
Target	HPV16
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03129
Expiration Date	12 months from date of receipt.