

Product Datasheet

FIPI (orb181391)

Description	Potent phospholipase D (PLD) inhibitor (IC50 values are 20 and 25 nM for PLD2 and PLD1 respectively).
Note	For research use only
Formula	$C_{23}H_{24}FN_5O_2$
MW	421.4764
CAS Number	[939055-18-2]
SMILES	<chem>C1=CC=CC2NC(=O)N(C3CCN(CCNC(=O)C4NC5C=CC(F)=CC=5C=4)CC3)C1=2</chem>
Expiration Date	12 months from date of receipt.